

TNOA – COVID UPDATE DRESS CODE FOR CLINICS

RECEPTION

- Full cover cloth dress
- N95 mask
- * Cloth mask over that
- *Face shield

FEVER OP STAFF

- Full cover cloth dress
- ❖ N95 mask
- Cloth mask over that
- Face shield
- Gloves if BP/Temperature are checked







NORMAL OP STAFF

- ❖ 3 ply/Cloth mask
- Face Shield
- Gloves if BP/Temperature Are checked



HOUSE KEEPING

- Plastic apron
- Elbow level gloves
- Leggings
- Cloth mask
- Face shield







COMMON

- ❖ Daily all their clothes, cloth mask to be soaked in 0.5% Sodium Hypochlorite solution for 15minutes.
- Then immediately to be washed with Soap & Water.
- Dry in open terrace.
- Face shield to be dipped in Soap water for 10 minutes and to be washed& dried to be used next day
- N-95 mask to be hanged in the room to be used after 72 hours or



keep the mask in hot air oven at 65 degree centigrade for 30 minutes remove and can be used.

Note: Preparation of 0.5% Sodium Hypochlorite solution 1Scoop Bleaching Powder in 10 Litres of Water.

Dr.V.ThangaveluPresident TNOA

Dr.V.Madhavan Hony. Secretary IMA FSS